

Quantum-yield high and low with large spot size →

High and low powers measured 2 spot sizes using lenses and flip mirror plus beam splitter to attenuate the light.

Author: Jean Matias

UCNP Sample: s013

spot sizes: 500um, 100um

Dye Sample: s017

Acquired on: 2021/04/30

Analysed on: 2021/05/04 (May the fourth be with you)

Last Updated: 2022/03/01

- Metadata available at the end of the notebook,
- The final data that generates the plots and their metadata are saved at the end of the notebook in `../data/analysed-data/` folder.

In [1]:

```
%matplotlib notebook
import sys
sys.path.append('../scripts/')
from lab import Sample, Analysis, BeamProfile
from calibration import *
```

For details on the analysis protocol and calculation:

```
>>> help(Analysis)
>>> help(Sample)
```

In [21]:

```
help(Analysis)
|   1. Get steady data - it averages data points after the equipment
and samples stabilised.
|   2. Calibrate the power meter output (pm) - Volts to mW using standard curves for that.
|       This is the transmitted power through the cuvette.
|   3. Calibrate the DAQ output (trigger) - Volts to mW using empty measurements.
|       This is the laser power, i.e. the power before reaching the cuvette.
|   4. Get the power density in the center - It uses the power before and after the cuvette
|       along the beam width.
|   5. Get absorbance - Sample-Absorbance is the absorbance of the whole sample.
|       On the other hand Absorbance is the absorbance of the compound without its diluter.
|       It's necessary to have the diluter absorbance in order to calculate this property.
|   6. Get absorbed power - Using absorbance it calculates the total power absorbed by the
```

In [3]:

```
help(Sample)
| Methods defined here:

|     __init__(self, path, expId, sampleType)
|         Initialize self. See help(type(self)) for accurate signature.

|     absorbance(self, diluterAbsorbance=None, filterQuery='laser_power >
| 0', recalculate=False)
|         Calculates absorbance of the emitter compound present on the aqueo
us medium.

|     Args:
|         diluterAbsorbance (float): Diluter absorbance used as a refere
nce to calculate the emitter's absorbance
|         recalculate (bool): If True it will recalculate the column

|     Returns:
|         pandas.DataFrame, float, float, str: Data, Absorbance mean, st
d and success, failure message.

|     absorbedpower(self, region=[0.45, 0.55], useMean=True, filterQuery='la
```

Selecting samples:

1. ucnp
2. diluters: e.g. water
3. empties
4. dyes
5. diluters: e.g. ethanol
6. empties

In [2]:

```
on20210504 = Analysis.datalog['exp_id'].str.contains('20210331') & ~Analysis.datalog['saved_name'].str.contains('sample=(?:s013|empty|water)/.*976nm')
Analysis.datalog[on20210504]
```

Out[2]:

	exp_id	saved_name	out_ch	range_start	range_end	range_step_size	step_reset
645	20210331-090934	./data/raw-data/qy_20210331_090934.csv	ao0=Laser	0.0	5.4	0.08	False
646	20210331-091102	./data/raw-data/qy_20210331_091102.csv	ao0=Laser	0.0	5.4	0.08	False
647	20210331-091239	./data/raw-data/qy_20210331_091239.csv	ao0=Laser	0.0	5.4	0.08	False
648	20210331-091659	./data/raw-data/qy_20210331_091659.csv	ao0=Laser	0.0	5.4	0.08	False
650	20210331-	./data/raw-	ao0=l aser	0.0	5.4	0.08	False

Beam Profile analysis

In [3]:

```
bp2 = BeamProfile('s2_20210430_976nm_140mA_ND1_ET03.png')
bp2_01 = BeamProfile('s2_20210430_976nm_180mA_ND1_BS50_ET03.png')


bp1 = BeamProfile('s1_20210430_976nm_100mA_ND1_ET002.png')
bp1_01 = BeamProfile('s1_20210430_976nm_100mA_ND1_BS50_ET003.png')

bp1.trim(4.5)
bp1_01.trim(4.5)
bp2_01.trim(1.5)
bp2.trim(1.2)

bp1.removebackground(2)
bp1_01.removebackground(2)
bp2_01.removebackground(3)
bp2.removebackground(3)
```

In [4]:


```
bp1.plotsurf()
```



In [11]:

```
bp1_01.plotsurf()
```


```
<IPython.core.display.Javascript object>
```



In [6]:

```
bp2_01.plotsurf()
```


```
<IPython.core.display.Javascript object>
```



In [8]:

```
bp2.plotsurf()
```

<IPython.core.display.Javascript object>



In [4]:


```
print(f'BP1: {bp1.beamwidth():.5f}, BP1_01: {bp1_01.beamwidth():.5f}, BP2: {bp2.beamwidth()}
```

BP1: 0.01113, BP1_01: 0.01060, BP2: 0.05300, BP2_01: 0.05221 in cm

In [15]:

```
_ = bp1.histprofile(show=True)
```


```
<IPython.core.display.Javascript object>
```



In [73]:

```
_ = bp2_01.histprofile(show=True)
```


```
<IPython.core.display.Javascript object>
```



In [74]:

```
_ = bp2.histprofile(show=True)
```

```
<IPython.core.display.Javascript object>
```



UCNP data

In [5]:

```
uData = Analysis()
```

In [6]:

```
# Beam S. size = s1 ~100um
uData.loaddata('ucnp', dataIds='20210430-120627', diluterId='20210430-120923', emptyId='202
# Beam S. size = s1 ~100um BS50:50
uData.loaddata('ucnp', dataIds='20210430-120134', diluterId='20210430-120020', emptyId='202

# Beam S. size = 400um
uData.loaddata('ucnp', dataIds='20210430-113345', diluterId='20210430-113227', emptyId='202
# Beam S. size = 400um BS50:50
uData.loaddata('ucnp', dataIds='20210430-114054', diluterId='20210430-114329', emptyId='202
```

In [7]:

```
uData.get('setbeamprofile', dict(beamprofileObj=bp1_01), which=[0, 1, 2])
uData.get('setbeamprofile', dict(beamprofileObj=bp1_01), which=[3,4,5])

uData.get('setbeamprofile', dict(beamprofileObj=bp2), which=[6, 7, 8])
uData.get('setbeamprofile', dict(beamprofileObj=bp2_01), which=[9,10,11])
```


20210430-120828: empty - Success!
20210430-120923: diluter - Success!
20210430-120627: ucnp - Success!
20210430-115026: empty - Success!
20210430-120020: diluter - Success!
20210430-120134: ucnp - Success!
20210430-112928: empty - Success!
20210430-113227: diluter - Success!
20210430-113345: ucnp - Success!
20210430-114220: empty - Success!
20210430-114329: diluter - Success!
20210430-114054: ucnp - Success!

In [10]:

3	115026	empty	empty	min	instead ND	BS50	976nm	600mW
4	20210430-120020	diluter	water	min	BS50:45 instead ND	BS50	976nm	600mW
5	20210430-120134	ucnp	s013	min	BS50:45 instead ND	BS50	976nm	600mW
6	20210430-112928	empty	empty	min	realigned	0	976nm	600mW
7	20210430-113227	diluter	water	min	realigned	0	976nm	600mW
8	20210430-113345	ucnp	s013	min	realigned	0	976nm	600mW
9	20210430-114220	empty	empty	min	BS50:45 instead ND	BS50	976nm	600mW
10	20210430-114329	diluter	water	min	BS50:45 instead ND	BS50	976nm	600mW
11	20210430-	210	.	.	BS50:45	BS50	976	600 mW

In [8]:

```
fig, axs = uData.view(includeDiluter=False, includeEmpty=False, label='spot_size')
```



Get Steady Data


In [9]:

```
uData.get('steadydata', dict(fromPt=400, toPt=1000, precision=2), includeDiluter=True, incl
```

```
20210430-120828: empty - Success!
20210430-120923: diluter - Success!
20210430-120627: ucnp - Success!
20210430-115026: empty - Success!
20210430-120020: diluter - Success!
20210430-120134: ucnp - Success!
20210430-112928: empty - Success!
20210430-113227: diluter - Success!
20210430-113345: ucnp - Success!
20210430-114220: empty - Success!
20210430-114329: diluter - Success!
20210430-114054: ucnp - Success!
```

In [10]:

```
fig, axs = uData.view(x='trigger', yList=['pm', 'apd'], includeDiluter=True, includeEmpty=T
```



Removing background

In [11]:

```
uData.get('removebackground', dict(channel='apd', baseline='trigger < 0.2'), includeEmpty=T
```

```
20210430-120828: empty - Success!
20210430-120923: diluter - Success!
20210430-120627: ucnp - Success!
20210430-115026: empty - Success!
20210430-120020: diluter - Success!
20210430-120134: ucnp - Success!
20210430-112928: empty - Success!
20210430-113227: diluter - Success!
20210430-113345: ucnp - Success!
20210430-114220: empty - Success!
20210430-114329: diluter - Success!
20210430-114054: ucnp - Success!
```

In [12]:


```
uData.get('removebackground', dict(channel='pm', baseline='trigger < 0.2'), includeEmpty=True)
```

20210430-120828: empty - Success!
20210430-120923: diluter - Success!
20210430-120627: ucnp - Success!
20210430-115026: empty - Success!
20210430-120020: diluter - Success!
20210430-120134: ucnp - Success!
20210430-112928: empty - Success!
20210430-113227: diluter - Success!
20210430-113345: ucnp - Success!
20210430-114220: empty - Success!
20210430-114329: diluter - Success!
20210430-114054: ucnp - Success!

In [13]:

```
fig, axs = uData.view(x='trigger', yList=['pm', 'apd'], includeDiluter=True, includeEmpty=True)
```


Figure 1



Power meter calibration

In [14]:

```
coeffs, r_val = calibratePM(laser=976, pmRange=600, degree=1)
```



In [15]:

```
coeffs
```

Out[15]:

```
[305.16892805223625, -0.37612596104816115]
```

In [17]:

```
pmCalibfunc600 = lambda x: coeffs[0] * x + coeffs[1]
```


In [18]:

```
uData.get('calibrate', dict(calibFunction=pmCalibfunc600, channel='pm', recalculate=True),
```

```
20210430-120828: empty - Success!
20210430-120923: diluter - Success!
20210430-120627: ucnp - Success!
20210430-115026: empty - Success!
20210430-120020: diluter - Success!
20210430-120134: ucnp - Success!
20210430-112928: empty - Success!
20210430-113227: diluter - Success!
20210430-113345: ucnp - Success!
20210430-114220: empty - Success!
20210430-114329: diluter - Success!
20210430-114054: ucnp - Success!
```

In [19]:


```
fig, axs = uData.view(x='trigger', yList=['transmitted_power', 'apd'], label='spot_size')
```



Trigger calibration

In [20]:

```
coeffs, r_val = uData.sample(0).polyfit(x='trigger', y='transmitted_power', degree=12, cond  
triggerCalib = lambda x: np.poly1d(coeffs)(x)
```



In [21]:


```
uData.get('calibrate', dict(calibFunction=triggerCalib, channel='trigger', recalculate=True)
```

20210430-120828: empty - Success!
20210430-120923: diluter - Success!
20210430-120627: ucnp - Success!

In [22]:

```
coeffs, r_val = uData.sample(3).polyfit(x='trigger', y='transmitted_power', degree=12, cond  
triggerCalib = lambda x: np.poly1d(coeffs)(x)
```

$$1.68e+02 x^5 + -2.03e+02 x^4 + 1.46e+02 x^3 + -48.2 x^2 + 32.2 x^1 + -11.1$$
$$R^2 = 0.999972$$



In [23]:


```
uData.get('calibrate', dict(calibFunction=triggerCalib, channel='trigger', recalculate=True)
```

```
20210430-115026: empty - Success!  
20210430-120020: diluter - Success!  
20210430-120134: ucnp - Success!
```

In [24]:

```
coeffs, r_val = uData.sample(6).polyfit(x='trigger', y='transmitted_power', degree=12, cond  
triggerCalib = lambda x: np.poly1d(coeffs)(x)
```

$$+02 x^5 + -4.64e+02 x^4 + 3.44e+02 x^3 + -1.33e+02 x^2 + 67.9 x^1 + -19.7$$
$$R^2 = 0.999974$$




In [25]:

```
uData.get('calibrate', dict(calibFunction=triggerCalib, channel='trigger', recalculate=True)
```

20210430-112928: empty - Success!
20210430-113227: diluter - Success!
20210430-113345: ucnp - Success!

In [26]:

```
coeffs, r_val = uData.sample(9).polyfit(x='trigger', y='transmitted_power', degree=12, cond=True)
triggerCalib = lambda x: np.poly1d(coeffs)(x)
```




In [27]:

```
uData.get('calibrate', dict(calibFunction=triggerCalib, channel='trigger', recalculate=True))
```

20210430-114220: empty - Success!
20210430-114329: diluter - Success!
20210430-114054: ucnp - Success!

In [28]:

```
fig, axs = uData.view(x='laser_power', yList=['transmitted_power', 'apd'], label='sample',
```



Power density at the centre of the cuvette

In [24]:

```
help(Sample.powdensatcentre)
```

Help on function powdensatcentre in module lab:

```
powdensatcentre(self, beamWidth=None, recalculate=False)
    Calculates power density at the centre of the cuvette.
    Power at the centre/beam area.
```

Args:

 recalculate (bool): if True it recalculates the power at the centre
of the cuvette.

In [177]:

```
help(Sample.poweratcentre)
```

Help on function poweratcentre in module lab:

```
poweratcentre(self, recalculate=False)
    Calculates power at the centre of the cuvette.
    Given by $\sqrt{P_o * P}$
```

Args:

```
    recalculate (bool): if True it recalculates the power at the centre
    of the cuvette.
```


In [29]:

```
uData.get('powdensatcentre', dict(recalculate=True), includeDiluter=True, includeEmpty=True)
```

```
20210430-120828: empty - Success!
20210430-120923: diluter - Success!
20210430-120627: ucnp - Success!
20210430-115026: empty - Success!
20210430-120020: diluter - Success!
20210430-120134: ucnp - Success!
20210430-112928: empty - Success!
20210430-113227: diluter - Success!
20210430-113345: ucnp - Success!
20210430-114220: empty - Success!
20210430-114329: diluter - Success!
20210430-114054: ucnp - Success!
```

In [30]:

```
fig, axs = uData.view(x='power_at_centre', yList=['transmitted_power', 'apd'], label='spot_
```



APD calibration

In [31]:

```
# Coefficient obtained from Apd-calibration notebook
coeffs = [0.00525246, 0] # dye reabsorption included
apdCalibfuncMin = lambda x: np.poly1d(coeffs)(x)
```

In [32]:

```
uData.get('calibrate', dict(calibFunction=apdCalibfuncMin, channel='apd', recalculate=True))
```

```
20210430-120627: ucnp - Success!
20210430-120134: ucnp - Success!
20210430-113345: ucnp - Success!
20210430-114054: ucnp - Success!
```


In [33]:

```
uData.get('emittedpower', dict(scatteringAtEmission=0.5, recalculate=True))
```

```
20210430-120627: ucnp - Success!
20210430-120134: ucnp - Success!
20210430-113345: ucnp - Success!
20210430-114054: ucnp - Success!
```

In [34]:

```
fig, axs = uData.view(x='pow_dens_at_centre', yList=['transmitted_power', 'emitted_power'],
```



Absorbance + Scattering

Sampleabsorbance represents the absorption coefficient plus scattering.
Scattering is obtained with white light measurements.


In [35]:

```
uData.get('sampleabsorbance', dict(recalculate=True), includeDiluter=True)
```

```
20210430-120923: diluter - Success!
20210430-120627: ucnp - Success!
20210430-120020: diluter - Success!
20210430-120134: ucnp - Success!
20210430-113227: diluter - Success!
20210430-113345: ucnp - Success!
20210430-114329: diluter - Success!
20210430-114054: ucnp - Success!
```

In [36]:

```
fig, axs = uData.view(x='time_daq', yList=['sample_absorbance'], label='sample', includeDil
```



In [37]:

```
for k in [1,4,7,10]:
    _, abMean, abStd, _ = uData.sample(k).sampleabsorbance(filterQuery='time_daq > 10', rec
print(f'{uData.sample(k)._dataID}-{uData.sample(k)._sampleType}: {abMean=}, {abStd=}'')
```

```
20210430-120923-diluter: abMean=0.5802187177682696, abStd=0.0039005442047185
99
20210430-120020-diluter: abMean=0.5750551746784278, abStd=0.0032505589472063
574
20210430-113227-diluter: abMean=0.5655113426351803, abStd=0.0013235354316370
31
20210430-114329-diluter: abMean=0.5702465562633373, abStd=0.0010001993915492
037
```

In [38]:

```
uData.get('absorbance', dict(recalculate=True), includeDiluter=True)
```

```
20210430-120923: diluter - Success!
20210430-120627: ucnp - Success!
20210430-120020: diluter - Success!
20210430-120134: ucnp - Success!
20210430-113227: diluter - Success!
20210430-113345: ucnp - Success!
20210430-114329: diluter - Success!
20210430-114054: ucnp - Success!
```

In [31]:


```
uData.sample(11).absorbance(filterQuery='time_daq > 10', recalculate=True)
```

Out[31]:

```
(0      NaN
 1      NaN
 2      NaN
 3      NaN
 4      NaN
 ...
63    0.343571
64    0.344981
65    0.346035
66    0.345985
67    0.344327
Name: absorbance, Length: 68, dtype: float64,
0.34002209804944966,
0.002499419719701178,
'20210430-114054: ucnp - Success!')
```

In [39]:

```
fig, axs = uData.view(x='laser_power', yList=['absorbance'], label='sample', includeDiluter
```




In [40]:

```
uData.get('absorbedpower', dict(filterQuery='time_daq > 10', recalculate=True, absorptionCo
```

```
20210430-120627: ucnp - Success!
20210430-120134: ucnp - Success!
20210430-113345: ucnp - Success!
20210430-114054: ucnp - Success!
```

In [41]:

```
fig, axs = uData.view(x='power_at_centre', yList=['absorbed_power', 'emitted_power'], label
```



Get kappa κ


In [42]:

```
uData.get('kappa', dict(period=10, recalculate=True))
```

```
20210430-120627: ucnp - Success!
20210430-120134: ucnp - Success!
20210430-113345: ucnp - Success!
20210430-114054: ucnp - Success!
```

In [43]:

```
fig, axs = uData.view(x='pow_dens_at_centre', yList=['emitted_power', 'kappa'], label='spot')
```



Quantum Yield ϕ


In [44]:

```
uData.get('quantumyield', dict(recalculate=True))
```

```
20210430-120627: ucnp - Success!
20210430-120134: ucnp - Success!
20210430-113345: ucnp - Success!
20210430-114054: ucnp - Success!
```


In [45]:

```
fig, axs = uData.view(x='absorbed_power', yList=['emitted_power', 'quantum_yield'], which=[
```



In [46]:

```
fig, axs = uData.view(x='pow_dens_at_centre', yList=['absorbed_power', 'emitted_power', 'qu
```



In [40]:

```
pl.tight_layout()
```

Fitting emission according to a simplified model

$$E = 2\phi_b \alpha_{976} A_{px} \frac{\lambda_e}{\lambda_a} \sum_k^N \frac{\rho_k^2}{\rho_b + \rho_k}$$


Being α_{976} the absorption coefficient of the UCNPs at 976 nm, A_{px} the area of a pixel, ρ power density

In [47]:

```
def formatcoefs(popt, pcov):
    perr = np.sqrt(np.diag(pcov))
    phib = popt[0]; phib_err = perr[0]
    rhob = popt[1]; rhob_err = perr[1]
    print(f'phi_b = {popt[0]:.4f}({perr[0]:.4f}), rho_b = {popt[1]/1000:.2f}({perr[1]/1000:.2f})')
    return(phib, phib_err, rhob, rhob_err)
```

In [67]:

```
popt, pcov, emissionF8, _ = uData.sample(8).fitemission(bounds=[[0, 0], [1, 1e+5]), show=True  
phib, phib_err, rhob, rhob_err = formatcoefs(popt, pcov)
```



```
phi_b = 0.0189(0.0004), rho_b = 38.22(1.50)W/cm^2
```

Converting units

In [69]:

```
rhob /= 1000  
rhob_err /= 1000  
phib *= 100  
phib_err *= 100
```

In [77]:

```
print(f"{{rhob:.1f} pm {rhob_err:.1f}, {phib:.2f} pm {phib_err:.2f}}")
```

```
rhob=38.2 pm rhob_err=1.5, phib=1.89 pm 0.04
```

Fitting graphs

In [81]:

```
pl.rcParams.update({'font.size': 13})
pl.rcParams.update({'font.family':'serif'})
def qy(p, phib, rhob): return (phib * p / (rhob + p))

fig = pl.figure(figsize=(7,8.5))
xlabel = 'Power density (W/cm$^2$)'
ylabel = 'Quantum Yield ($\$\%)'


ax1 = pl.subplot2grid((3,2), (0,0))
ax2 = pl.subplot2grid((3,2), (0,1))
ax3 = pl.subplot2grid((3,2), (1,0))
ax4 = pl.subplot2grid((3,2), (1,1))
ax5 = pl.subplot2grid((3,2), (2,0), colspan=2)

ax1.set_xlabel(xlabel)
ax2.set_xlabel(xlabel)
ax3.set_xlabel(xlabel)
ax4.set_xlabel(xlabel)
ax5.set_xlabel(xlabel)
ax1.set_ylabel(ylabel)
ax2.set_ylabel(ylabel)
ax3.set_ylabel(ylabel)
ax4.set_ylabel(ylabel)
ax5.set_ylabel(ylabel)

# fig, ax = pl.subplots()
ax4.plot(uData.sample(2).data()['pow_dens_at_centre'], uData.sample(2).data()['quantum_yiel'])
ax3.plot(uData.sample(5).data()['pow_dens_at_centre'], uData.sample(5).data()['quantum_yiel'])
ax2.plot(uData.sample(8).data()['pow_dens_at_centre'], uData.sample(8).data()['quantum_yiel'])
ax1.plot(uData.sample(11).data()['pow_dens_at_centre'], uData.sample(11).data()['quantum_yiel'])

ax5.plot(uData.sample(2).data()['pow_dens_at_centre'], uData.sample(2).data()['quantum_yiel'])
ax5.plot(uData.sample(5).data()['pow_dens_at_centre'], uData.sample(5).data()['quantum_yiel'])
ax5.plot(uData.sample(8).data()['pow_dens_at_centre'], uData.sample(8).data()['quantum_yiel'])
ax5.plot(uData.sample(11).data()['pow_dens_at_centre'], uData.sample(11).data()['quantum_yiel'])

# fitted model
power = np.arange(0, 1932, 0.5)
ax5.plot(power, qy(power, 2*phib, rhob), '-.', color='red', label='Beam compensation')
pl.ylabel(ylabel)
pl.xlabel(xlabel)
pl.legend(fontsize=11)
fig.tight_layout()
```



In [82]:

```
# save figure
fig.savefig('../plots/Quantum_Yield_pannel.png')
```

In [58]:

```
# emission
data = '../data/emission/s013_20210521_emission_QEP025961.txt'
data = pd.read_csv(data, sep='\t', skiprows=15, names=['Wavelength', 'Intensity'])
```

In [59]:

```
data
```

Out[59]:


	Wavelength	Intensity
0	346.573	1.25
1	347.359	-0.75
2	348.146	-1.75
3	348.933	246.25
4	349.719	227.25
...
1038	1123.340	397.25
1039	1124.047	-3.75
1040	1124.754	-5.75
1041	1125.460	-2.75
1042	1126.167	-0.75

1043 rows × 2 columns

In [60]:

```
pl.rcParams.update({'font.size': 15})
fig, ax = pl.subplots(figsize=(7, 6))
ax.plot(data['Wavelength'], data['Intensity']/data['Intensity'].max(), color='k', label='Emission')

ax.set_xlabel('Wavelength ($nm$)')
ax.set_ylabel('Emission intesity ($a.u.$)')
ax.set_ylim([0,1.05])
ax.set_xlim([400,900])
ax.yaxis.tick_right()
ax.yaxis.set_label_position("right")
# ax.grid()
pl.tight_layout()
```



In [50]:

```
fig.savefig('../plots/spectrum.pdf')
```

Exporting analysed data that generate the plots

In [83]:

```
for idx in uData.details().index:  
    fileName = '../data/analysed-data/' + uData.details().loc[idx, 'exp_id'] + '_' + uData.  
    metadataName = '../data/analysed-data/' + uData.details().loc[idx, 'exp_id'] + '_meta  
  
    uData.sample(idx).data().to_csv(fileName)  
    uData.sample(idx)._metadata.to_csv(metadataName)
```

In [84]:

```
uData.sample(11)._metadata
```

Out[84]:

	field_name	data_type	data_format	example	standard_units	plot_label	des
0	time	float	$^{([0-9]?[.][1-9][\d]*.) ([0-9]*$)}$	123423.982734	seconds	Time (s)	nu
1	trigger	float	$^{[-]?([0-9]?[.] [1-9][\d]*.) ([0-9]*$)}$	-1.998340	volts	Trigger (V)	C
2	pm	float	$^{[-]?([0-9]?[.] [1-9][\d]*.) ([0-9]*$)}$	-1.998340	volts	PM (V)	C
3	apd	float	$^{[-]?([0-9]?[.] [1-9][\d]*.) ([0-9]*$)}$	-1.998340	volts	APD (V)	C
4	time_daq	float	$^{[-]?([0-9]?[.] [1-9][\d]*.) ([0-9]*$)}$	0.029380	s	Time (s)	Dac
							Tra
5	transmitted_power	float	$^{[0-9]^*[.][0-9]^*}$	5.500000	mW	Power (mW)	(ca
6	laser_power	float	$^{[0-9]^*[.][0-9]^*}$	5.500000	mW	Power (mW)	(fro ca
7	power_at_centre	float	$^{[0-9]^*[.][0-9]^*}$	100.565000	mW	Power (mW)	tr
8	pow_dens_at_centre	float	$^{[0-9]^*[.][0-9]^*}$	10000.565000	W/cm ²	Power (W/cm ²)	d tr
							Ap
9	apd_power	float	$^{[0-9]^*[.][0-9]^*}$	5.500000	mW	Power (mW)	- (ca
10	emitted_power	float	$^{[0-9]^*[.][0-9]^*}$	10.565000	mW	Emission (mW)	poi di
							Abs
11	sample_absorbance	float	$^{0.[0-9]^*}$	0.400000	(cm ⁻¹)	Absorbance (cm ⁻¹)	obt log wl

	field_name	data_type	data_format	example	standard_units	plot_label	des
12	absorbance	float	^0.[0-9]*	0.400000	cm^{-1}	Absorbance (cm^{-1})	Abs obt abs
13	absorbed_power	float	^0.[0-9]*	0.400000	mW	Abs. pow. (mW)	C
14	kappa	float	^[0-9]*.[0-9]*	10000.565000	-	κ	E slo loc
15	quantum_yield	float	^[0-9]*.[0-9]*	10000.565000	-	ϕ	C y lam lai
16	power_density	float	^[0-9]*.[0-9]*	5.500000	W/cm^2	Power (W/cm^2)	C frc
17	absorption	float	^[0-9]*.[0-9]*	5.500000	mW/cm^2	Absorption (mW/cm^3)	Abs A p vol
18	emission	float	^[0-9]*.[0-9]*	5.500000	mW/cm^2	Emission (mW/cm^3)	E p vol

In []: